

To customers using MeGUARD.

Thank you very much for using MeGUARD.

Mounting Position of the MeGUARD Frame Please pay attention to the position where the frame is placed.	
Correct Wearing Example	Incorrect Wearing Example
	

Before use, please read the following instructions carefully and use the product correctly.

This product is designed to be properly worn at the usual position where a headband is placed (on the top of the head) to prevent it from easily falling off.

If worn in an incorrect position (such as the forehead), the frame may fall off during use.

The frame itself is not designed for permanent use. If the frame has expanded in shape or lost its elasticity due to deterioration compared to a new frame, we recommend replacing the frame.

